






# Hemosep<sup>®</sup>

## Product Catalogue

	Product	Product Code	Quantity per Box
	<b>Hemosep Shaker</b>	<b>CR4213</b>	<b>1</b>
	<b>Hemosep Adult Pack</b> Hemosep Adult 1000ml new membrane cell concentrator & blood collection bag.	<b>CR4426</b>	<b>1 set</b>
	<b>Hemosep Paediatric Pack</b> Hemosep Paediatric 500ml new membrane cell concentrator & blood collection bag.	<b>CR4427</b>	<b>1 set</b>
	<b>Hemosep Intra Operative Kit: CONTAINS;</b> CR4421: Hemosep Blood Collection Reservoir (manufactured and packed/ sterile by SIDD India) CR4434: Hemosep Suction Kit	<b>CR4428</b>	<b>1 set</b>
	<b>Hemosep Suction Kit: CONTAINS;</b> One A/A Line with attached suction wand and including slider On/ Off valve One vacuum connection line One Spare Suction wand (Tubing)	<b>CR4434</b>	<b>1 set</b>
	<b>Hemosep One to Three Connector</b>	<b>CR4432</b>	<b>1</b>
	<b>Hemosep Intra-Operative Stand</b>	<b>CR4429</b>	<b>1</b>